

## TYPICAL PRODUCT PARAMETERS

	<b>Invacare Transport (TRANI7FR, TRANI9FR)</b>	<b>Lightweight Transport (LTTB17FR, LTTRI9FR, LTTI319FR)</b>	<b>Bariatric Transport (TRHD22FR)</b>	<b>Lite Transport (ALRI9HBFR, ALBI9HBFR)</b>
<b>SEAT WIDTH (BETWEEN ARMS)</b>	14 <sup>1</sup> / <sub>2</sub> and 16 <sup>1</sup> / <sub>2</sub> inches	14 <sup>1</sup> / <sub>2</sub> and 16 <sup>1</sup> / <sub>2</sub> inches	22 inches	17 inches
<b>SEAT DEPTH</b>	16 inches	16 inches	18 inches	16 inches
<b>OVERALL WIDTH</b>	21 to 23 inches	21 to 23 inches	28 inches	23 <sup>1</sup> / <sub>2</sub> inches
<b>OVERALL DEPTH (WITH-OUT RIGGINGS)</b>	24 <sup>1</sup> / <sub>4</sub> inches	24 <sup>1</sup> / <sub>4</sub> inches	29 inches	28 inches
<b>SEAT-TO-FLOOR</b>	19 <sup>1</sup> / <sub>2</sub> inches	19 <sup>1</sup> / <sub>2</sub> inches	19 <sup>1</sup> / <sub>2</sub> inches	19 <sup>1</sup> / <sub>2</sub> inches
<b>BACK STYLE</b>	Fold-down	Fold-down	Fold-down	Fold-down
<b>BACK HEIGHT</b>	18 inches	16 inches	19 inches	17 inches
<b>ARM STYLES</b>	Permanent	Permanent	Removable	Permanent
<b>FRONT RIGGINGS</b>	Swingaway footrest	Swingaway footrest	Swingaway footrest	Swingaway footrest
<b>REAR AXLE</b>	Permanent	Permanent	Permanent	Permanent
<b>REAR WHEELS</b>	8 X 1 inch	8 X 1 inch	12 <sup>1</sup> / <sub>2</sub> X 2 <sup>1</sup> / <sub>4</sub> inch	12 X 1 <sup>3</sup> / <sub>4</sub> inch
<b>WHEEL LOCKS</b>	Push to lock	Push to lock	Push to lock	Push to lock
<b>HAND BRAKES</b>	N/A	N/A	Push down to lock, pull up to stop	Push down to lock, pull up to stop
<b>CASTER SIZE</b>	8 X 1 inch (Standard)	8 X 1 inch (Standard)	8 X 1 inch (Standard)	8 X 1 inch (Standard)
<b>UPHOLSTERY</b>	Black Nylon	Black Nylon	Black Nylon	Black Nylon
<b>WEIGHT (APPROXIMATE)</b>	26 lbs (without front riggings)	19 lbs (without front riggings)	41 lbs (without front riggings)	19 lbs (without front riggings)
<b>SHIPPING WEIGHT (APPROXIMATE)</b>	30 lbs (without front riggings)	26 lbs (without front riggings)	45 lbs (without front riggings)	26 lbs (without front riggings)
<b>WEIGHT LIMITATION</b>	250 lbs	250 lbs	400 lbs	250 lbs